

*Otto Bock*<sup>®</sup>

QUALITY FOR LIFE



# SiliconGel Liner 6Y75



with SKINGUARD<sup>®</sup> Technology

Product Information





# SiliconGel Liner 6Y75




## Skin-friendly – antibacterial – neutral odor

The 6Y75 SiliconGel Liner with SKINGUARD® Technology include the new and effective additive Sanitized®. The high-tech functionality of antibacterial additive Sanitized® protects your liner against bacterial involvement and thus also against bad odors. It preserves the material, improves hygiene and the odor reduction.



# SiliconGel Liner 6Y75 with SKINGUARD® Technology



(Sanitized® inhibits the growth of bacteria on polymer surfaces)

The natural antibacterial effect of silver against a broad-range of bacteria is undisputed and scientifically accepted. It damages the cell membranes of bacteria and inhibits their growth. In addition, it hinders the reproduction of the bacteria responsible for odor development.

Independent lab tests according to the JIS Z 2801 Standard have proven that 6Y75 SiliconGel Liners treated with Sanitized® caused a 95% reduction of Staphylococcus aureus (gram positive).

Like all Otto Bock SKINGUARD® Technology products, Sanitized® distinguishes itself through its dependability and brand strength. Sanitized® was subjected to internationally recognized risk assessments for people and the environment.





Technical modifications reserved.

**Otto Bock®**


**QUALITY FOR LIFE**

**Otto Bock HealthCare GmbH**

Max-Näder-Straße 15 · 37115 Duderstadt/Germany · Phone +49 5527 848-1424 · Fax 05527 72330 · [export@ottobock.de](mailto:export@ottobock.de) · [www.ottobock.com](http://www.ottobock.com)

© Otto Bock · 646A215=GB - 06.08

# Convincing – Features & Benefits



## + Antibacterial additive Sanitized®

To protect and improve your wellbeing, the 6Y75 SiliconGel Liner contains an antibacterial additive that prevents bacterial growth and bad odors in your liner effectively and without irritating the skin.




## + Soft silicone, with a silky, skin-friendly internal surface.

Improved skin condition and comfort while not losing the adhesion for suspension.




## + Distal elongation reduction system that does not interfere with circumferential stretch

Our patented integrated finger matrix minimizes longitudinal movement (pistoning) without restricting the radial elongation (around the limb).



## + Ultra-smooth and durable silver outer textile

Easier donning and doffing of the liner especially for geriatrics and patients with lower hand involvement.



## + Soft distal cap

The soft distal cap is slightly flexible, thus making it easier to step into and engage with the prosthesis.



# SiliconGel Liner 6Y75 with SKINGUARD® Technology

The SiliconGel Liner with SKINGUARD® Technology are suitable for amputees that require a durable but soft silicone in the low to moderate activity levels and perform best in a specific weight bearing socket design.

## Winning Combinations

The 6Y75 SiliconGel liner is ideal for shuttle system sockets in combination with 6A30 or 6A20 shuttle locks.

The liner is available in sizes: 160, 180, 200, 210, 220, 235, 250, 265, 280, 300, 320, 340, 360, 380, 400.

